

# QIAgility pipetting robot

## What is it?

- A robot for molecular laboratory liquid handling.

## Where is it?

- Molecular Biology Lab, rm 220.1, Agr Ed Building (W077)

## What does it do?

- Automated PCR reaction set-up;
- Automated DNA, RNA and other solution serial dilutions;

## Who can use it?

- Anyone who has been inducted and is undertaking UNE related research, teaching or outreach preparation.

## What training and access?

- You will need to be inducted in use of the QIAgility and also for access to the lab that houses it. Access after hours is not permitted.

## Who to see?

- Contact technician: John Pesor;
- Contact academic: Shubiao Wu [swu3@une.edu.au](mailto:swu3@une.edu.au); 67732238

