

2100 Bioanalyzer

What is it?

- DNA, RNA, and protein sizing, quantitation and quality control system.

Where is it?

- Molecular Biology Lab, rm 220.1, Agr Ed Building (W077)

What does it do?

- Size DNA, RNA and protein;
- Quantify DNA, RNA and protein;
- Measure the integrity of DNA and RNA.

Who can use it?

- Anyone who has been inducted and is undertaking UNE related research, teaching or outreach preparation.

What training and access?

- You will need to be inducted in use of the Bioanalyzer and also for access to the lab that houses it. Access after hours is not permitted.

Who to see?

- Contact technician: John Pesor;
- Contact academic: Shubiao Wu swu3@une.edu.au; 67732238

